

MA 39, MA 40, MA 41 & MA 42

Versatile, Portable Audiometers



Select the model that best suits your needs - from basic air conduction to a combination of air-bone or air-bone-speech to a complete two-channel air-bone-speech audiometer with master hearing aid capabilities.

- Quiet, reliable electronic switches
- Insert phones capability
- Built-in RS232 interface for PC and NOAH compatibility
- Bright, easy to read display screen
- Lightweight, portable soft case

EXCELLENCE IN HEARING SCREENING

SINCE 1937



MA 39, MA 40, MA 41 & MA 42

Hard case



Soft case



MA 39

The MA 39 is a portable audiometer that performs quick, accurate pure tone threshold testing.

MA 40

The MA 40 provides pure tone threshold as well as air and bone testing. This unit is NOAH compatible and includes a built-in RS232 interface to transfer data to a computer.

MA 41

The MA 41 is a portable audiometer that provides air, bone and speech testing and has the capacity to drive two free field speakers without an additional amplifier. This unit is NOAH compatible and includes a built-in RS232 interface to transfer data to a computer.

MA 42

The MA 42 provides air, bone and speech testing as well as monaural and binaural master hearing aid testing which makes it ideal for hearing aid fitting and dispensing. It includes separate insert phone jacks, CD/tape input, monitor speaker and monitor headset. Powerful built-in amplifiers provide ample power to connect directly to your speaker system and save you the expense of buying external power amplifiers. This unit is NOAH compatible and includes a built-in RS232 interface to transfer data to a computer.

With the master hearing aid you can select four commonly used frequency slopes to approximate hearing aid acoustics and let your patients experience the difference hearing aids can make.

If air-bone-speech is air-bone-speech... why pay thousands more?

Technical data

Part	MA 39	MA 40	MA 41	MA 42
<i>Standard Accessories</i>				
TDH 39 headset	■	■	■	■
B-71 bone receiver and cord		■	■	■
Bone headband		■	■	■
Monitor/mic headset			■	■
Audiogram pad	■	■	■	■
Operation Manual	■	■	■	
Patient response switch				■
<i>Optional Accessories</i>				
Patient response switch	■	■	■	
Audiocup headset	■	■	■	■
Insert phones	■	■	■	■
Talkback mic			■	■
Sound room kit	■	■	■	■
Free field kit			■	■
Sound field kit			■	■
Patch cord	■	■	■	■

CPT Codes

	MA 39	MA 40	MA 41	MA 42
92551 Screening, pure tone	■	■	■	■
92552 Threshold, air	■	■	■	■
92553 Air and bone		■	■	■
92555 Threshold, speech			■	■
92556 Threshold/discrimination			■	■
92557 Basic comprehensive			■	■
92590 Hearing aid exam and selection, monaural				■
92591 Binaural				■
92592 Hearing aid check, monaural				■
92593 Binaural				■

Specifications subject to change without notice



EXCELLENCE IN HEARING SCREENING

SINCE 1937

