

MI 34

Middle Ear Analyzer



Maico combines the speed and simplicity of an impedance screener with the added functions of reflexes and reflex decay tests that are usually found only on more expensive diagnostic equipment.

- Ipsi and contra reflex testing
- Reflex decay testing
- Intact Eustachian tube function testing
- Perforated Eustachian tube function testing
- Optional 1 kHz high frequency probe tone
- Adjustable comfort probe
- LCD graphical display
- Quiet, fast internal printer

EXCELLENCE IN HEARING SCREENING

SINCE 1937



MI 34

Automatic test results, quickly and easily

Diagnostic functionality

The Maico MI 34 performs automatic impedance measurements for middle-ear evaluations. Besides the typical acoustic reflex tests used in hearing screening, the MI 34 provides diagnostic functions such as ipsi and contra reflex decay as well as intact and perforated Eustachian tube function. There is also the option to include a high frequency (1000 Hz) probe tone function in the instrument which is beneficial when testing neonates.

Quality features

The MI 34 was designed for flexibility and easy operation. Users can customize test protocols, default settings and printouts. It includes a large backlit LCD screen for enhanced views of testing in progress and an internal thermal printer that provides test results in seconds. In addition, the instrument's probe can be set at three

different angles to easily maneuver around the patient for a complete seal. The MI 34 is lightweight and can include an optional soft-sided carrying case for a portable office solution.



MI 34 with accessories

CPT codes

92567	Tympanometry
92568	Acoustic reflex
92569	Acoustic reflex decay

Specifications subject to change without notice

ONE YEAR LIMITED WARRANTY: This warranty is extended to the original purchaser of the instrument, by Maico, through the Distributor from whom it was purchased. The warranty covers defects in material and workmanship for a period of one year from date of delivery of the instrument to the original purchaser. Accessories which are purchased from Maico at the same time as the instrument are warranted for one year from the date of purchase. For additional information contact your Maico Distributor.

EXCELLENCE IN HEARING SCREENING

SINCE 1937

