

NEW FROM WELCH ALLYN...

## TWICE THE ON-TIME, HALF THE WEIGHT



### **INTRODUCING THE WELCH ALLYN POWER HANDLE ...SMARTER TECHNOLOGY**

For the first time in a power handle, Welch Allyn introduces smarter Lithium Ion technology to provide twice the on-time with just half the weight. The Lithium Ion Power Handle's new Lithium Ion technology powers any of our 3.5v instrument heads for 110 minutes – twice as long as the conventional handle. With an updated, ergonomic design and weighing only 4.4 ounces, the handle is more comfortable for any size hand. An indicator light warns you when the battery is running low, and your choice of charging modules makes recharging simple and convenient. Made of durable CYCOLOY® construction, the Lithium Ion Power Handle will provide years of reliable quality.

**WelchAllyn®**  
Lithium Ion  
Power Handle

# LITHIUM ION TECHNOLOGY PROVIDES TWICE THE ON-TIME WITH JUST HALF THE WEIGHT.

## Lithium Ion Power Handle Features

- Half the weight of conventional 3.5v handles
- Average 110 minutes of on-time between charging
- Smallest diameter of any full-size 3.5v power handle
- Built to last using CYCOLOY® resin
- Low-battery indicator
- Lithium Ion battery technology
- ETL listed to UL 2601 and CAN/CSA C22.2 No. 601.1 medical electrical equipment standards

## Lithium Ion Battery Technology

- No “memory” effect
- Low battery indicator
- No need to discharge between charging cycles



Soft grip rheostat intensity control section

Textured handle

Battery section durable CYCOLOY® blend construction

Low-battery indicator light (Handle shown at actual size)

## CHOICE OF CHARGING MODULES



### Removable, Pull-out AC Charging Module

- Easily removable from handle
- Retractable AC prongs
- Can be attached or detached to handle during exams



### Lithium Ion Handle Well Adapter

- Fits Welch Allyn 71110 desk charger
- “Right sizing” top adjusts to different well heights

\*CYCOLOY is a registered trademark of General Electric Company



*Charger module can be detached when not in use.*



*Charger module can be left attached for enhanced convenience and comfort for larger hands.*



*Retractable module can be plugged into any 110v outlet.*

## ORDERING INFORMATION

- 71900** Lithium Ion Power Handle with AC Charging Module
- 71930** Lithium Ion Power Handle with Well Adapter
- 71911** Handle without AC Charging Module
- 71920** Lithium Ion Well Adapter only
- 71950** AC Charging Module only
- 71960** Replacement battery

U.S. Patent No. 6,130,520

For more information about the Welch Allyn Lithium Ion Power Handle, contact Welch Allyn at **888-359-2394**.

**WelchAllyn®**  
Lithium Ion  
Power Handle

4341 State Street Road, Skaneateles Falls, NY 13153-0220 Tel: (888) 359-2394, (315) 685-2679 Fax: (315) 685-2174  
www.welchallyn.com infomedical@mail.welchallyn.com