

# The Newest & Best Compact WDVRA Video Wireless Pediatric Testing System

The Wireless WDVRA Visual Reinforcement Audiometry (VRA/COR) system provides endless "instant on" child-friendly digital animations to use either alone or in conjunction with built in audio stimulus. System includes HD 10" widescreen LCD with "instant on" technology.

**Quick  
& Easy to  
Install**

**Best  
Value!**

**Small  
&  
Compact**

## Key Features

- Bright LCD Technology Display
- 4 GB CF Memory with Continuous Cycling
- Preprogrammed Wireless Remote Control
- "Instant On" Technology
- Tilt & Swivel Mounting Brackets

## Optional Accessories

- Centering Monitor
- Assistant Wireless Remote
- Rolling Stand
- Upgradeable CF Memory
- Customized Systems Available

**Bright  
LCD  
Display!**

**All Inclusive  
Package!**



Call for more information

**NORTHEASTERN TECHNOLOGIES GROUP**

Leaders in Hearing & Balance Testing Equipment

**800-229-4634**



**AVAILABLE NATIONWIDE!**



# ***The Wireless WDVRA Video Visual Reinforcement Audiometry (VRA/COR) System***

## **Specifications**

- Screen size : 10.2" (Diagonal)
- Active area : 222.7(H) x 130.5(V)
- Resolution : 720 dpi HD
- Brightness : 220 cd/m
- Contrast ratio : 400 : 1
- Power input : 12V DC, 2A
- Power consumption : 24W
- Remote controller : 9V DC

## **Warranty**

The WDVRA is warranted for a period of one year from the date of purchase against all manufacturing defects and workmanship. Defective products may be returned to the factory with pre authorization for repair or replacement at the manufacturer's discretion. Any improper use, tampering or negligent use of product is also subject to discretion of the manufacturer as a warranty claim. Batteries are not covered by warranty.

## **WDVRA-20**

Includes: One WDVRA-20 wireless controller, two 10 inch 720 dpi widescreen monitors each containing a formatted 4GB CF. Two adjustable swivel & tilt arm wall mount bracket for flat panel. Two 12vdc 2A power adapters with extensions.

## **WDVRA-30**

Includes: One WDVRA-30 wireless controller, three 10 inch 720 dpi widescreen monitors each containing a formatted 4GB CF. Three adjustable swivel & tilt arm wall mount bracket for flat panel. Three 12vdc 2A power adapters with extensions. (Third monitor used for centering stimuli).

## **NORTHEASTERN TECHNOLOGIES GROUP**

*Leaders in Hearing & Balance Testing Equipment*

40 Glen Street, Glen Cove, New York 11542

**800-229-4634**

Fax: 516-674-3507

**www.northeasterntech.com**

